

QAQC UD DOSER

Accurate Unit Dose Sampling

Ideal for accurately sampling specific volumes of liquids. Designed specifically for unit dose sampling of liquids in the pharmaceutical industry.

- Suitable for sampling low to medium viscosity liquids.
- Ideal for low volume (unit dose) sampling.
- Made from Sampler body - 316 stainless steel,
Syringe - PP & PE

Sample Chamber is designed so that the outside of the syringe is not contaminated during sampling. The syringe can therefore be used as a carrying container for taking the sample back to the laboratory.

UD DOSER SOP

Operation



1. Select a new syringe



2. Connect the syringe to the push rod



3. Attach outer cover. The sampler is now ready for use.



4. Pull up handle to take the sample, the graduations allow an accurate sample to be taken

5 ml Doser

675 5010A600 600mm
675 5010A1200 1200mm
675 5010A1800 1800mm

Spare 5 ml Syringes

675 5010S05 Pack of 100

60 ml Doser

675 5010C600 600mm
675 5010C1200 1200mm
675 5010C1800 1800mm

Spare 60 ml Syringes

675 5010S60 Pack of 25

QUICK QUOTE

ONLINE PRICING

