


Part No.	055H1000S
Description	Disposable ViscoThief – 1000mm
Material of Construction:	
Tip	HDPE (virgin, FDA approved – CFR 177.1520)
Outer Tube	HDPE (virgin, FDA approved – CFR 177.1520)
Inner Rod	HDPE (virgin, FDA approved – CFR 177.1520)
Clutch	HDPE (virgin, FDA approved – CFR 177.1520)
Piston Housing	HDPE (virgin, FDA approved – CFR 177.1520)
Piston Seal	Silicone Rubber (virgin, FDA approved – CFR 177.2600)
Handle	HDPE (virgin, FDA approved – CFR 177.1520)



Method of Construction:	Ultrasonic welding
Nominal Length	
Fully closed:	991mm
Fully extended:	1818mm
Diameter of Outer Tube:	21mm
Diameter of hole in tip:	10mm
Nominal Weight of Sampler:	158g
Assembly & Packing Environment:	Class 100,000 Medical Cleanroom
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Samplers per Box:	20
BSE/TSE Free:	Yes

Recommended Storage Conditions: Dry and ambient temperature